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ENGLISH

COMPOSITION:

Geistlich Mucograft® / Geistlich Mucograft® Seal is a fully resorbable collagen matrix obtained by standardized, controlled manufacturing processes. The matrix is made of porcine collagen without further cross-linking. The collagen is extracted from veterinary certified pigs and is carefully purified. Geistlich Mucograft® / Geistlich Mucograft® Seal is sterilized in double packaging by Gamma-irradiation.

Geistlich Mucograft® / Geistlich Mucograft® Seal has two structures:

The compact structure consists of compact collagen which allows tissue adherence as a prerequisite for favorable wound healing. This layer has a smooth texture with appropriate elastic properties to accommodate suturing. A second structure consists of a thick, porous collagen spongy scaffold. This spongy surface is placed next to the host tissue to allow tissue integration.

PROPERTIES / ACTIONS:

The low antigenicity and excellent biocompatibility favor the use of Geistlich Mucograft® / Geistlich Mucograft® Seal in dental surgery. The longfibred microstructure of Geistlich Mucograft® / Geistlich Mucograft® Seal readily absorbs fluid. The coherent collagen fibers of Geistlich Mucograft® / Geistlich Mucograft® Seal form a basic tissue matrix. As a result, the matrix adheres well to the surrounding and underlying tissues. Inflammatory reactions have not been observed, but cannot be excluded. Geistlich Mucograft® / Geistlich Mucograft® Seal is approximately 2.5–5.0 mm thick. Fixation by sutures or pins is possible.

INDICATIONS:

Geistlich Mucograft® and / Geistlich Mucograft® Seal are indicated for:

- > covering of implants placed in immediate or delayed extraction sockets
- > localized gingival augmentation to increase keratinized tissue (KT) around teeth and implants

- > alveolar ridge reconstruction for prosthetic treatment
- > recession defects for root coverage.

INSTRUCTIONS FOR USE:

The general principles of sterile handling and patient medication must be followed when using Geistlich Mucograft® / Geistlich Mucograft® Seal.

- > The defect is exposed by means of (a) properly prepared flap(s) and the usual surgical procedures in creating a properly prepared site are undertaken.
- > When using Geistlich Mucograft® / Geistlich Mucograft® Seal for extraction socket management, the defect is filled with a space maintaining material, such as autologous bone or bone substitute (e.g. Geistlich Bio-Oss Collagen®). Such defects should not be overfilled.
- > The Geistlich Mucograft® / Geistlich Mucograft® Seal matrix is cut, if needed, to the required size and shape with surgical scissors. Geistlich Mucograft® / Geistlich Mucograft® Seal must remain dry during trimming and application to the defect.
- > The compact structure must face outwards, away from underlying bone, and the spongy structure must face toward the bone. In the presentation of 8 mm diameter, the spongy structure is striped for easier differentiation of the two sides. The striped spongy structure must face the bone.
- > The matrix is applied over the prepared site and stabilized using sutures or the flap. Excessive pressure should be avoided as it may compress the matrix.
- > Complete penetration of the matrix by blood and exudates allows close adaptation and adhesion of the matrix to the underlying surface.
- > Due to the compact structure of this matrix, fixation or suturing is possible. Fixation or suturing of the matrix may be indicated, depending on the nature of the particular defect and to avoid displacement of the matrix. After placement in covered sites, the mucoperiosteal flap should be sutured over the matrix without tension.

SPECIAL INSTRUCTIONS FOR USE IN PERIODONTOLOGY:

A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planing, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success.

POST-OPERATIVE CARE:

The usual postoperative care and medication should be given. Further prosthetic treatment should only be carried out after a healing period to ensure complete soft tissue regeneration. Previous studies have indicated that soft tissue healing, i.e., wound closure and the resolution of normal inflammation, occurs within 4–8 weeks. But clinical judgment, taking into account the healing of both soft and hard tissues in the treated patient, should be applied before prosthetic treatment. In case of bacterial contamination rinsing with appropriate bactericidal solutions is recommended.

In the rare event that matrix removal is necessary, the tissues adjacent to the matrix should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent

to the residual matrix. Following careful reflection of the surrounding tissue, the remaining portion of the matrix can be excised and the area debrided to remove any inflamed or infected tissue.

LIMITATIONS FOR USE:

Contraindications:

Geistlich Mucograft® / Geistlich Mucograft® Seal should not be placed where symptomatic infection or inflammation exists. Before placement, the surgeon should be confident that any symptomatic or recent infection has been properly treated. Patients with known allergy to collagen should not be treated with Geistlich Mucograft® / Geistlich Mucograft® Seal.

Precautions:

Geistlich Mucograft® / Geistlich Mucograft® Seal should be used with special caution in patients with autoimmune diseases or uncontrolled metabolic diseases (e.g. diabetes, osteomalacia, thyroid disorder), as well as in case of a prolonged corticosteroid therapy or radiotherapy in the oral cavity. In addition it is not recommended to use the matrix in more than one layer. Healing of soft tissues might be compromised in patients with insufficiently vascularized soft tissues (e.g. in smokers). The material has not been tested in pregnant or lactating women or in children.

Side effects:

In addition the following potential side effects, resulting from the surgical process, may arise after placing a matrix into the oral cavity: soft tissue dehiscence, hematoma, increased sensitivity and pain, redness and inflammation.

STORAGE AND HANDLING:

Do not use after the expiration date. The content of the double packaging is designed for **single use only**. Do not reuse or resterilize Geistlich Mucograft® / Geistlich Mucograft® Seal. The matrix should be handled using sterile gloves or sterile instruments. The matrix is sterile unless the package has been opened, damaged, or otherwise contaminated. Store in a dry place at room temperature (15–25 °C/ 59–77 °F).

PRESENTATION:

Geistlich Mucograft® / Geistlich Mucograft® Seal is packed in sterile double packaging. Each double packaging contains one matrix.

Re-Order no.	Sizes
50-1520	Geistlich Mucograft® 15 x 20 mm
50-2030	Geistlich Mucograft® 20 x 30 mm
500024	Geistlich Mucograft® 30 x 40 mm
50-0808	Geistlich Mucograft® Seal 8 mm diameter
500726	Geistlich Mucograft® Seal 12 mm diameter

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG
Bahnhofstrasse 40
CH-6110 Wolhusen
Switzerland
Made in Switzerland

This product is sold and distributed under US Patent No. 6,713,085

Caution:

Federal law restricts this device to sale by or on the order of a dentist or physician.

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