Geistlich Bio-Gide® Perio

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:
Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

de GERMAN

en ENGLISH

es SPANISH

fr FRENCH

Geistlich Bio-Gide® Perio

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Rev. 2021-12

Geistlich Bio-Gide® Perio

Resorbable Bilayer Collagen Membrane for Guided Tissue and Bone Regeneration

Geistlich

Biomaterials

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:
Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

ENGLISH

Geistlich Bio-Gide® Perio is a biocompatible, sterile resorbable bilayer collagen membrane for guided tissue and bone regeneration.

COMPOSITION:

Geistlich Bio-Gide® Perio is a pure collagen membrane obtained by standardized, controlled manufacturing processes. The collagen is extracted from veterinary certified pigs and is carefully purified to avoid antigenic reactions. Geistlich Bio-Gide® Perio is sterilized in double blisters by γ -irradiation. Geistlich Bio-Gide® Perio has a bilayer structure: The porous surface — facing the bone — will allow the ingrowth of the bone forming cells. The dense surface — facing the soft tissue — will prevent the ingrowth of fibrous tissue into the bone defect. The membrane is made of collagen without further cross-linking.

PROPERTIES / ACTION:

The low antigenicity and its biocompatibility allows the use of Geistlich Bio-Gide® Perio in dental surgery. Animal studies have shown that the collagen is resorbed in bone cavities within 24 weeks. The fibrous microstructure of Geistlich Bio-Gide® Perio is hydrophilic. The membrane retains its structural integrity even when wet. Fixation by sutures or pins is possible. Geistlich Bio-Gide® Perio has a smooth dense surface. As this prolongs the water uptake, cutting the membrane for periodontal defects is facilitated. The coherent collagen fibers swell and form a unified basic tissue structure. As a result, the adaption to the bone wall and the complete closure of the bony defects is easily achieved. Inflammatory reactions have not been observed, but due to the special physical properties and elongated resorption time a certain inflammatory reaction after placing the membrane cannot be excluded.

INDICATIONS FOR USE:

Geistlich Bio-Gide® Perio is intended for the following uses:

- > augmentation around implants placed in immediate extraction sockets
- > augmentation around implants placed in delayed extraction sockets

- > localized ridge augmentation for later implantation
- > alveolar ridge reconstruction for prosthetic treatment
- > filling of bone defects after root resection, cystectomy, removal of retained teeth
- > guided bone regeneration in dehiscence defects
- > guided tissue regeneration procedures in periodontal defects.

CONTRAINDICATIONS:

Geistlich Bio-Gide® Perio should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated. Geistlich Bio-Gide® Perio should not be used in patients with a known allergy to collagen.

DIRECTIONS FOR USE:

The general principles of sterile handling and patient medication must be followed when using Geistlich Bio-Gide® Perio:

- > The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures (e.g., curettage) are undertaken.
- > The bone defect is filled with a space-making material, e.g., autogenous bone, Geistlich Bio-Oss®, Geistlich Bio-Oss Collagen®, etc. The defect should not be overfilled.
- > The size of Geistlich Bio-Gide® Perio is adapted by scissors according to the exposed defect. The Geistlich Bio-Gide® Perio package contains sterile stencil templates in 4 different basic periodontal shapes. The membrane should overlap the walls of the defect by at least two millimeters to allow complete bone contact and to prevent gingival connective tissue invasion below the material.
- > The smooth surface should face the soft tissue, and the rough side should face the bone.
- > Geistlich Bio-Gide® Perio is applied over the defect without further treatment and held in place with moderate pressure. The period of time necessary to apply pressure will vary with the degree of bleeding. Adhesion to the bone surface is facilitated by the hydrophilic properties of the membrane.
- > Complete penetration of the membrane by blood and exudate allows adaptation and adhesion of the membrane to bony structures.
- > Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination.
- > Due to the high tensile strength of this membrane, fixation is possible. Fixation may be indicated to avoid membrane displacement due to loading or mobilization.
- > The mucoperiosteal flap can be sutured over the collagen membrane (e.g., single sutures and deep mattress sutures).
- > Complete wound closure is not required for Geistlich Bio-Gide® Perio when used in extraction sites and alveolar ridge defects with an average vertical defect length and alveolar ridge width up to 5.7 mm and 18 mm, respectively. In all other cases, complete wound closure should be attempted.

Please note: Although it was not observed for Geistlich Bio-Gide® Perio, it cannot be totally excluded that in case of membrane exposure the ridge width gain may be limited. However, ridge width gain noted for early membrane exposures still was sufficient to allow for the placement of dental implants.

- > In order to avoid the formation of a long junctional epithelium when treating periodontal defects, the coronal portion of Geistlich Bio-Gide® Perio should be adapted closely to the tooth (e.g. with additional fixation using suture material).

POST-OPERATIVE CARE:

- > The patient should be monitored closely.
- > If the membrane becomes exposed (in case of closed healing), the dehiscence usually heals by itself within several weeks.
- > Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions is recommended.
- > In the event that early membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.
- > To allow for undisturbed bone regeneration underneath the membrane surgical reentry should not be done before 4 to 6 months postoperatively.

PRECAUTIONS:

Geistlich Bio-Gide® Perio should only be used by trained dentists or oral surgeons.

In case of membrane exposure during the healing phase the resorption time may be accelerated. A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success.

Do not implant the stencil templates. Their perforations enable them to be distinguished from the membrane.

If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability. The long-term safety and effectiveness of the membrane in conjunction with bone-filling materials has not yet been established.

Absolute stability of the membrane is important for guided bone regeneration and is a vital condition for therapeutic success, and the smallest movement on the tissue underneath is to be avoided.

Effect on pregnant or lactating women is not known.

Geistlich Bio-Gide® Perio should be used with special caution in patients with any medications or clinical conditions that negatively influence tissue healing especially:

- > acute or chronic infection at the surgical site
- > uncontrolled metabolic diseases (e.g., diabetes, osteomalacia, thyroid disorder)
- > prolonged corticosteroid therapy
- > autoimmune diseases
- > radiotherapy
- > heavy smoking.

Local use of the membrane with medicinal products, alcohol, disinfectants or antibiotics has not been studied.

ADVERSE REACTIONS:

As Geistlich Bio-Gide® Perio is a collagen product in very rare cases allergic reaction may occur. Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, and local inflammation.

STORAGE AND HANDLING:

If the sterile packaging (double sterile barrier system) is damaged or opened, the product must not be used. The contents of the double blister is designed for single use only.

Do not resterilize Geistlich Bio-Gide® Perio. Resterilization is not possible. Discard any unused material. Do not use after expiry date.

Store at controlled room temperature 15°–25° C (59°–77° F) and in a dry place. Keep away from sunlight.

The membrane should be handled using sterile gloves or sterile instruments.

PRESENTATION: Geistlich Bio-Gide® Perio is packed in sterile double blisters.

Re-order #	Size
20148	1 Geistlich Bio-Gide® Perio membrane, 16 x 22 mm

Geistlich Perio-System Combi-Pack

Re-order #	Size
20147	1 Geistlich Bio-Gide® Perio membrane, 16 x 22 mm 1 Geistlich Bio-Oss Collagen®, 100 mg (1 block 100 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

 **Manufacturer:**

Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

Made in Switzerland

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.

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